

INNOVATOR INSIGHT

Scaling-up the production of stem cell-derived extracellular vesicles in stirred-tank bioreactors

Aurélie Tacheny

Exosomes are extracellular vesicles (EVs), that are secreted by different cells of the body and are considered as important players in cell-to-cell communication. Exosomes, and more specifically stem cell-derived exosomes are currently of great interest as cell-free therapeutic tools due to their strong diagnostic and therapeutical potential in various diseases models. To explore their use in various biomedical frameworks, large amounts of high-quality exosomes need to be produced. To meet this challenge, cell culture in the controlled environment offered by stirred-tank bioreactor can facilitate the standardized expansion of large numbers of viable cells and enhance the reproducible production of extracellular vesicles. In this article, we discuss the challenges of exosome production and provide insights in its scale-up in stirred-tank bioreactors.

Cell & Gene Therapy Insights 2021; 7(8), 1077–1083

DOI: 10.18609/cgti.2021.140

THE POTENTIAL OF EXOSOME TECHNOLOGY IN ADVANCED THERAPIES

EVs can be broadly classified, based on their mechanism of release and size, in exosomes

(less than 150 nm in diameter), microvesicles/microparticles/ectosomes (100–1000 nm) and apoptotic bodies (>1000 nm) [1]. Originally considered as cellular waste product, exosomes are a population of naturally occurring

mobile extracellular vesicles released by different cell types *in vivo* under various physiological and pathological conditions. They are capable of transporting various bioactive molecules such as proteins, lipids and different types of nucleic acids. Exosomes play a role in mediating intracellular communication, in addition to modulating immunoregulatory processes, tumor metabolism, regenerative and degenerative processes. Beyond their use as biomarkers, exosomes offer promising opportunities for therapeutic applications, as vehicles for drug delivery, but also as cell-free therapeutic tools [2,3].

THE CHALLENGES OF EXOSOME MASS PRODUCTION

The preclinical and clinical development of exosome technology as a cell-free therapeutic and drug delivery platform requires large quantities of exosomes. The therapeutic doses required will vary depending on the specific purpose of the study and methods by which the exosomes are quantified. In earlier preclinical trials, 10^9 to 10^{11} particles were administered per animal to achieve a biological outcome [4]. Dosing is highly variable and will depend on the therapeutic target and method of delivery [5].

Mass production of exosomes is challenging for several reasons. A delicate balance must be achieved in terms of quantity and quality. The final product must be pure and retain its biological properties. To achieve this, manufacturing processes must be scalable, controlled and standardized to improve reproducibility. Unfortunately, most of the current methods for enrichment, isolation and characterization of exosomes do not meet these criteria.

The exosome production workflow is a multi-step process composed of three main steps:

1. The upstream process, in which the cells used as the exosome source are expanded and exosomes are enriched in the medium;

2. The downstream process, during which the secreted exosomes are separated from the other components of the conditioned medium; and finally
3. The quality control phase which is based on exosome characterization.

Considering the first step of the exosome production workflow, the upstream process leading to efficient exosome secretion, it is necessary to highlight some key factors, in order to ensure a successful and reproducible upstream process.

Extracellular vesicles are varied [6] and cell culture parameters, like the cell seeding density, passage number, cell confluence, and medium composition, can have an impact on the exosomes secreted quantity and quality [7-9]. To achieve robust proliferation, culture media are usually supplemented with fetal bovine serum (FBS) as a source of growth factors. However, these animal-derived supplements contain large amounts of exogenous exosomes that contaminate the vesicles of interest. To avoid such contamination, xeno-free cell culture media can be used during the cell expansion and the vesicle collection steps.

The growth support on which the cells are expanded is also important. Conventional 2D culture carriers for adherent cells, such as flasks or multilayer flasks, have a limited growth area per volume. Three dimensional, microcarrier-based cell culture in a bioreactor is more physiologically representative and easier to scale up. Several teams have already published data demonstrating that 3D culture is beneficial for exosome yield and activity [4,10]. The 3D culture parameters need to be optimized to ensure robust attachment to the microcarriers and subsequent cell growth. For example, different types of carriers are available, which differ among other things in the size, surface charge, and surface coating, and one needs to optimize the type of carriers to be used. The initial cell seeding volume, the ratio between the number of cells and microcarriers, and also the initial adhesion phase, for example in a reduced volume with

or without agitation need to be optimized as well. Culture parameters that may affect exosome production and quality include, pH, oxygen tension, temperature, agitation, and feeding strategy. All should be kept as stable as possible and the control of dissolved oxygen and temperature should aim at mimicking the physiological conditions to avoid cellular stress.

EXOSOME PRODUCTION IN STIRRED-TANK BIOREACTORS

Cultivation in bioreactors opens up new possibilities for process monitoring and control compared with conventional cell culture flasks. pH, temperature, and DO, among other parameters, can be monitored and controlled in real time, making it easier, on the one hand, to achieve conditions resembling the physiological situation, and on the other hand, to reproduce cultivation conditions from batch to batch.

The use of stirred-tank bioreactors offers several advantages, if an increase in the cell number is required. When using cell culture plates or flasks, increasing the number of cells usually involves 'scale-out', meaning increasing the number of culture vessels. Stirred-tank bioreactors also allow for 'scale-up', meaning increasing the culture volume by using larger


vessels. In a scale-up approach less vessels need to be handled compared to the scale-out situation, which can reduce manual work and save space. Stirred-tank design allows similar vessel geometries and capabilities at different scales, facilitating a smooth transition to larger volumes. It is the bioreactor design for which most research on scale-up phenomena has been conducted and it is well established in industrial cell culture bioprocessing.

With the following example we summarize, what kind of results can be expected when using stirred-tank bioreactors. For a detailed description of the study the reader is referred to a previous publication [11]. In this case, the DASbox® Mini Bioreactor System was used for the development of an exosome production process and the optimization of culture conditions at small scale, before a potential scale-up. Figure 1 shows a summary of the experimental strategy for the upstream process.

Mesenchymal stem cells derived from adipose tissue were selected as the exosome source. After thawing, the cells recovered for several days in flasks and the appropriate number was seeded on microcarriers in BioBLU® 0.3c Single-Use Bioreactors. Cell growth was monitored for 10 days. Growth factors were added 3 days after seeding. During the expansion phase, regular checks of cell adhesion and viability were


► FIGURE 1

Schematic representation of cell expansion and exosome production in BioBLU 0.3c Single-Use Bioreactors controlled with a DASbox Mini Bioreactor System.



► FIGURE 2

Isolation and characterization of stem cell-derived extracellular vesicles from hADSC-conditioned media.



Total protein profile of fractionated conditioned media (Nanodrop analysis) and CD63 determination by ELISA detection of SEC fractions.

performed. After 5 days of cell growth, the culture medium was completely replaced with collection medium. Exosome enrichment occurred during 2 additional days prior to harvest and subsequent downstream processing. All culture media used in this study were xeno-free [11].

The main process parameters used during this experiment have been previously summarized [11]. This paper includes information on the type of microcarriers, the inoculation cell density, the final working volume, the different physicochemical parameters used, and agitation.

The analysis of cell viability using fluorescent calcein-AM and fluorescence


microscopic evaluation revealed that during the five days post seeding cells grew and progressively formed small aggregates. Automated cell counting after trypsinization from the carriers confirmed the efficient proliferation with a very short lag phase. The final cell density was ca. 8.3×10^5 cells/mL, corresponding to $2.0^8 \times 10^8$ cells per bioreactor and a 34.7-fold increase in the cell number [11].

For downstream processing, the exosome conditioned medium was collected, centrifuged, filtered, and concentrated using Centriplus® YM 10 filter units (Amicon®). A small sample was taken for particle size distribution analysis. The remaining volume was used for size exclusion chromatography. More than 20 fractions were collected and the total amount of protein evaluated. Moreover, for each fraction we analyzed the presence of CD63 by ELISA and Western Blot analysis. As a negative control, a similar volume of non-conditioned medium has followed the same purification steps.

Figure 2 shows the results of the ELISA test. The total amount of protein in each fraction is indicated by the green line; the bars indicate the expression of an exosome marker (CD63). A high relative content of CD63-positive vesicles was found in the

► FIGURE 3

Particle size distribution of isolated hADSC-derived extracellular vesicles determined by dynamic light scattering analysis.



fractions 16 to 20. The peak in fraction 19 corresponds to a relative quantity of 9.1×10^{10} extracellular vesicles in this fraction. In the non-conditioned medium, an analogous pattern of total protein content was observed but no expression of the CD63 was detected. The presence of CD63 was confirmed by Western Blot analysis [11].

The particle size distribution was evaluated by dynamic light scattering analysis (DLA) (Figure 3). We found a mono-modal distribution of approximately 60 to 110 nm with a mean size of approximately 80 nm, which was the expected size for small extracellular vesicles [6].

CONCLUSION

In comparison with conventional 2D culture platforms, bioreactors offer important advantages to optimize the upstream process of an exosome production workflow among which process monitoring, control and scalability opportunities. Following the cell expansion phase, they ensure a fine control of the main bioprocess parameters during additional incubation days in EV collection medium, ensuring optimal conditions for EV quality preservation. This study is an example highlighting the potential of hADSCs and hADSC-derived extracellular vesicle production in a stirred-tank bioprocess system.

REFERENCES

1. Yanez-Mo M, Siljander PRM, Andreu Z *et al.* Biological properties of extracellular vesicles and their physiological functions. *J. Extracell. Vesicles* 2014; 4(1): 27066.
2. Hu Q, Su H, Li J *et al.* Clinical applications of exosome membrane proteins. *Precis. Clin. Med.* 2020; 3(1): 54–66.
3. Gurunathan S, Kang MH, Jeyaraj M, Qasim M, Kim JH. Review of the isolation, Characterization, biological function and multifarious therapeutic approaches of exosomes. *Cells* 2019; 8(4): 307.
4. Haraszti RA, Miller R, Stoppato M *et al.* Exosomes produced from 3D cultures of MSCs by tangential flow filtration show higher yield and improved activity. *Mol. Ther.* 2018; 26(12): 2838–47.
5. Shekari F, Nzari A, Kashani SA, Hajizadeh-Saffar E, Lim R, Baharvand H. Pre-clinical investigation of mesenchymal stromal cell-derived extracellular vesicles: a systematic review. *Cytotherapy* 2021; 23(4): 277–84.
6. Tkach M, Kowal J, Théry C. Why the need and how to approach the functional diversity of extracellular vesicles. *Phil. Trans. R. Soc.* 2017; 373: 20160479.
7. Patel D, Gray KM, Santharam Y, Laminchane TN, Stroka KM, Jay SM. Impact of cell culture parameters on production and vascularization bioactivity of mesenchymal stem cell-derived extracellular vesicles. *Bioeng. Transl. Med.* 2017; 2(2): 170–9.
8. Palviainen M, Saari H, Kärkkäinen O *et al.* Metabolic signature of extracellular vesicles depends on the cell culture conditions. *J. Extracell. Vesicles* 2019; 8(1): 1596669.
9. Gudbergsson JM, Johnsen KB, Skov MN, Duroux M. Systematic review of factors influencing extracellular vesicle yield from cell cultures. *Cytotechnology* 2016; 68(4): 579–92.
10. Thippabhotla S, Zhong C, He M. 3D cell culture stimulates the secretion of in vivo like extracellular vesicles. *Nature Sci. Rep.* 2019; 9(1): 13012.
11. Tejerina S, Dufey V, Hoet JF, Tachenay A, De Longueville F. The DASbox® Mini Bioreactor System as a tool for process development and stem cell-derived exosome production in standardized culture conditions. Eppendorf Application Note 440. 2021.

AFFILIATION

Aurélien Tachenay

Eppendorf Application Technologies
SA, Rue du Séminaire 20A, B-5000
Namur, Belgium

eppendorf

AUTHORSHIP & CONFLICT OF INTEREST

Contributions: All named authors take responsibility for the integrity of the work as a whole, and have given their approval for this version to be published.

Acknowledgements: None.

Disclosure and potential conflicts of interest: Aurélie Tacheny is employed by Eppendorf Application Technologies SA.

Funding declaration: The author received no financial support for the research, authorship and/or publication of this article.

ARTICLE & COPYRIGHT INFORMATION

Copyright: Published by Cell and Gene Therapy Insights under Creative Commons License Deed CC BY NC ND 4.0 which allows anyone to copy, distribute, and transmit the article provided it is properly attributed in the manner specified below. No commercial use without permission.

Attribution: Copyright © 2021 Eppendorf AG. Published by Cell and Gene Therapy Insights under Creative Commons License Deed CC BY NC ND 4.0.

Article source: Invited; externally peer reviewed.

Submitted for peer review: Aug 9 2021; **Revised manuscript received:** Sep 15 2021; **Publication date:** Sep 30 2021.

NEW: Bioprocess
Control System



Scale Up!

Seamlessly transition your process to the next level

Many engineering parameters need to be considered to ensure a successful tech transfer from bench to pilot and production scale. Efficient and reliable scale-up is no hurdle anymore with the BioFlo® 720 bioreactor control system*.

- > **Scalability:** Compatible with the Thermo Scientific™ HyPerforma™ 5:1 Single-Use Bioreactors (S.U.B.)
- > **Efficiency:** Intuitive software tools such as *Auto Calibrate* and *Auto Inflate*
- > **Application Driven:** Documentation packages available to help assist qualification into cGMP biomanufacturing environments



www.eppendorf.com/bioflo720

* BioFlo 720 is not a medical device as defined by the Food and Drug Administration or other regulatory authorities. Eppendorf®, the Eppendorf Brand Design, and BioBLU® are registered trademarks of Eppendorf AG, Germany. BioFlo® is a registered trademark of Eppendorf, Inc., USA. Thermo Scientific™ and HyPerforma™ are trademarks of Thermo Fisher Scientific Inc., USA. All rights reserved, including graphics and images. Copyright © 2021 by Eppendorf AG.